
Micro-Stain™ Bacterial Live/Dead Staining Kit (Fluorometric)

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- ▶ Fluorescence microscope / Fluorometric microplate reader / Flow Cytometry
- ▶ Vortexer / Centrifuge
- ▶ Tissue culture-treated 96well microwell plate, Dark clear-bottomed
- ▶ Pipette & Sterile tips / Microtube
- ▶ Slide glass / cover glass
- ▶ PBS (or 0.85% NaCl)
- ▶ Pure Isopropyl alcohol

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
Figure 1. SYBR Green/PI staining analyzed with fluorescence microscopy. Fluorescence images of the same samples at 530 nm (green, Live cells) for SYBR Green signal, 630 nm (red, Dead cells) for PI signal and merged images are shown. 50:50 ratio of live and dead cells of *E. coli* was used.


Figure 2. SYBR Green/PI staining can assess the viability of *E. coli*. The graph shows the percentage of live *E. coli* cells (G/R ratio), calculated based on fluorescence intensity ratios from SYBR Green (G; live) and PI (R; dead) signals.

Related products

- | | |
|-------------|--|
| - QB-0100 | Quanti-Micro™ Microbial Viability Assay Kit (Colorimetric) |
| - MS-BC1000 | Micro-Stain™ Bacterial CTC Staining Kit (Fluorometric) |
| - MS-BM1000 | Micro-Stain™ Bacterial Live/Dead Metabolism Assay Kit (Fluorometric) |
| - MS-YC1000 | Micro-Stain™ Yeast Viability Assay Kit (Fluorometric) |

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